

GENEWIZ Sanger Seq Drop Box Service for Moffitt Cancer Center!

Moffitt Cancer Center, Tampa, FL

- **Fast Turnaround:** Results posted to your account the Same Day samples are received by GENEWIZ
- **High Quality Data:** ABI3730xl Sanger Sequencers produce read lengths 750-1000 bp with Quality Score and Contiguous Read Length
- **Convenient:** User-friendly LIMS system for ordering, tracking, and retrieving results available; free shipping via GENEWIZ Drop Boxes
- **Superior Service:** Expert Technical Support available at 8 AM – 8 PM Eastern Time, Monday – Friday.
- **Cost Savings!:** Please call (346) 387-5613, or email Deidra.Conner@GENEWIZ.com for detailed pricing information



GENEWIZ offers three service levels for Sanger sequencing of purified plasmid and PCR templates: Pre-mixed, Pre-defined, and Custom.

If you need help choosing your optimal Sanger service level, please contact us:

[Email](#) | Phone (1-877-GENEWIZ (436-3949), Ext. 2) | [Live Chat](#).

How To Order

1. Visit <https://www.genewiz.com/> to **Register** an account or **Login** to your existing account. (banner, top of page)
2. Place your order online using your GENEWIZ account. Select "Sanger Sequencing", choose a template, and complete order details. Print the order receipt upon submission. You will need a PO or credit card to complete your order.
3. Prepare samples according to GENEWIZ [Sample Submission Guidelines](#) under "Resources".
4. No charge shipping - Convenient Sample Drop off and delivery via GENEWIZ Drop Box at:

MRC 1302

Latest Drop off Time: 5 PM
5. View your results from your online GENEWIZ account

For more information, please contact:

Deidra Conner
Sales Executive II
Cell: (346) 387-5613
Deidra.Conner@genewiz.com
www.genewiz.com

GENEWIZ, Inc.
115 Corporate Boulevard
South Plainfield, NJ 07080
877-GENEWIZ (436-3949)